

11.

AV Access Tunneler

PAGE 1

LOGO: Bard Peripheral Vascular

HEADLINE: AV ACCESS SHEATH TUNNELER SET

SUBHEAD: Advancing the Art of Tunneling and Graft Placement
with Accuracy and Interchangeability

PIC A: CU Hand Reaching in Tray for Tunneler
CAPTION: INITIATE

PIC B: CU Hand Assembling Tunneler
CAPTION: CREATE

PIC C: CU Hand Directing Tunneler
CAPTION: NAVIGATE

PIC D: CU Hand Placing Graft through Tunneler
CAPTION: FACILITATE

INSET PIC: Tunneler Tray Set

PAGE 2

LOGO: Bard Peripheral Vascular

HEADLINE: AV ACCESS SHEATH TUNNELER SET

SUBHEAD: Advancing the Art of Tunneling and Graft Placement
with Accuracy and Interchangeability

COPY: **Initiate Planning and Select Components**

- Tunneler set is intended to aid creating subcutaneous tunnels for placement of AV grafts up to 7 MM
- Choose from a selection of curved or straight components in 5 Or 6/7 mm diameters

Create Tunneler Assembly

- Easily attach the interchangeable handle on either end of the threaded inner rods
- When heeded, use multiple sheaths to precisely match anatomical needs

Navigate Tunneler with Accuracy and Interchangeability

- Tunnels can be formed from any direction
- Tunneler and suture attachment mechanism provide for numerous graft options

Facilitate Graft Placement with Confidence and Convenience

- Outer sheaths help to:
 - Prevent twisting or kinking.; protecting graft integrity
 - Confirm appropriate placement before removal
 - Permit free use of saline during graft insertion and reposioning

CALL: For ordering information, call 1-800-321-4254

EDITOR NOTE: NEED INPUT for CATALOG NUMBERS and DESCRIPTIONS

Order Info: TBD

Sign-Off: For more information, contact:
Bard Peripheral Vascular, Inc.
P.O. Box 11740
Tempe, AZ 85280
USA
TEL: 480-894- 9515
www.bardpv.com

Code numbers / reorder number